




KCARE SHOWER STOOL WITH ARMS - HEIGHT ADJUSTABLE, CLIP-ON SEAT



SKU: KA222ZA

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The KCare Shower Stool with Arms has a plastic seat that can easily be removed for cleaning/replacement and the height and angle can be adjusted. There are arms provide support to users while sitting or transferring. Rubber tips on each leg provide anti-slip stability. This shower stool is durable and is available with an epoxy powder coated zinc treated mild steel, white or aluminium frame and moulded PVC seat.

FEATURES

- Designed to assist those who have difficulty showering.
- Height adjustable legs .
- Plastic seat for ease of cleaning .
- Powder coated white 22mm steel tubing .

SPECIFICATIONS

Overall Width	510 mm
Overall Depth	380 mm
Product Weight	3 kg
Max User Weight	110 kg
Seat Width	340 mm
Seat Depth	290 mm
Min Seat Height	510 mm

Max Seat Height

630 mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au