



KCARE STANDARD UTILITY CHAIR - LOW BACK - HEIGHT AND WIDTH ADJUSTABLE



SKU: KA586KD-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

This Utility Chair from KCare features a contoured and padded backrest, padded armrests and height adjustable legs, to ensure the chair can be a comfortable seating solution for all.

Useful when a slightly taller chair is required to sit or stand up comfortably or after hip or knee surgeries, while the armrests also assist with pushing up from a seated position.

Can be used at the dining table, in the bedroom as a dressing chair, in the bathroom to sit on whilst getting dry or dressing/undressing or in the living room as a chair that will provide firm support at an easy height for activities or just to relax.

FEATURES

- The Utility chair is a multipurpose adjustable seat height chair
- Chair width and depth is adjustable with a total adjustment range of 60mm in both width and depth
- Padded arms and back
- Low back version has curved backrest
- Provides extra postural support
- Seat height is adjustable
- Rubber stoppers for stability
- Suitable for use in any area

SPECIFICATIONS

Min Overall Width	600 mm
Max Overall Width	660 mm
Overall Depth	510 mm
Min Seat Width	460 mm
Max Seat Width	510 mm
Seat Depth	460 mm
Min Seat Height	430 mm
Max Seat Height	580 mm
Backrest Height	440 mm
Max User Weight	200 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au