



KCARE STANDARD UTILITY CHAIR - LOW BACK - HEIGHT AND WIDTH ADJUSTABLE

~~\$357.00~~ Original price was: ~~\$357.00~~. ~~\$295.00~~ Current price is: **\$295.00**. Sale Available Online Only!



SKU: KA586KD-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

This Utility Chair from KCare features a contoured and padded backrest, padded armrests and height adjustable legs, to ensure the chair can be a comfortable seating solution for all.

Useful when a slightly taller chair is required to sit or stand up comfortably or after hip or knee surgeries, while the armrests also assist with pushing up from a seated position.

Can be used at the dining table, in the bedroom as a dressing chair, in the bathroom to sit on whilst getting dry or dressing/undressing or in the living room as a chair that will provide firm support at an easy height for activities or just to relax.

FEATURES

- The Utility chair is a multipurpose adjustable seat height chair
- Chair width and depth is adjustable with a total adjustment range of 60mm in both width and depth
- Padded arms and back
- Low back version has curved backrest
- Provides extra postural support
- Seat height is adjustable
- Rubber stoppers for stability
- Suitable for use in any area

SPECIFICATIONS

Min Overall Width	600 mm
Max Overall Width	660 mm
Overall Depth	510 mm
Min Seat Width	460 mm
Max Seat Width	510 mm
Seat Depth	460 mm
Min Seat Height	430 mm
Max Seat Height	580 mm
Backrest Height	440 mm
Max User Weight	200 kg