



## KCARE TILT IN SPACE SHOWER COMMODE WITH SWINGAWAY FOOTRESTS

**\$3,892.00**



**SKU: KA11507**

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The KCare Tilt In Space Shower Commode with Swingaway Footrests is a versatile and functional mobility aid designed for individuals with limited mobility. This commode features a tilt-in-space function that allows for smooth and safe transfers, as well as enhanced comfort during use. It also comes with swingaway footrests that can be adjusted to accommodate different leg lengths and provide added support during transfers. It comes with an extra high PVC strap backrest for improved comfort and safety, as well as padded armrests that can be swung up to allow for unrestricted side transfers. Additionally, the commode is equipped with 4x 125mm total locking castors that provide easy manoeuvrability, making it ideal for use in various settings.

**Please note: Padded seat is not included**

### FEATURES

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.
- Extra high PVC strap backrest for improved comfort and safety.
- 10-35° tilt angle.
- Integrated push-handle for manoeuvrability.
- Swing-up padded arms for unrestricted side transfer.
- Swingaway and removable height adjustable footrest with contoured padded calf supports.
- 4x 125mm total locking castors.
- Designed to be used with a K Care patented removable seats.
- **Padded seat not included** - [Open Front PU Seat](#), [Closed Front PU Seat](#), or [Handmade Open Seat](#) available.

- **[Pan](#), and [pan carrier](#) not included - sold separately**

## ACCESSORIES

- [Open Front PU Seat](#), [Closed Front PU Seat](#), or [Handmade Open Seat](#)
- [Commodes Pan and Bowl Carrier](#)
- [Commode Bowl with Lid and Handle](#)
- [Commode Bed Pan](#)

## SPECIFICATIONS

Overall Width	580 mm
Overall Depth	1070 mm
Overall Height	1453 mm
Seat Width	445 mm
Seat Depth	430 mm
Floor To Seat Height	540 mm
Armrest Type	Loop
Armrest Height	365 mm
Backrest Height	855 mm
Footrest Type	Swing Away
Product Weight	21 kg
Max User Weight	200 kg