



KCARE TRANSIT MOBILE SHOWER COMMODE WITH SLIDEOUT FOOTPLATE

\$1,105.00 - \$2,362.00Price range: \$1,105.00 through \$2,362.00



SKU: KA112

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

he Transit Mobile Shower Commode with Slideout Footplate is a versatile and practical solution for individuals with mobility difficulties. This commode features a slideout footplate that provides additional support for the user's feet and allows for easy side transfers. It is designed with a comfortable padded seat and a 3-strap backrest for optimal comfort and support. The commode is equipped with four 125mm total locking castors and an integrated push-handle for easy manoeuvrability. Additionally, it has swing-up arms that enable unrestricted side transfer and a removable pan for easy cleaning.

Please note: Padded seat is not included

FEATURES

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.
- 3 strap backrest for comfort and easy access.
- Integrated push-handle for manoeuvrability.
- Swing-up arms for unrestricted side transfer.
- Slide out/retractable non-weight bearing platform.
- 4x 125mm total locking castors.
- Designed to be used with a K Care patented removable seats.
- Available in 445, 500, 550, 600 and 650mm seat widths.
- Seat, pan, and pan carrier not included sold separately



SPECIFICATIONS

					Floor						Max
Seat	Overall	Overall	Overall	Seat	To	Armrest	Armrest	Backrest	Footrest	Product	User
Width	Width	Depth	Height	Depth	Seat	Type	Height	Height	Type	Weight	Weight
					Height						weight
445	480	840	1005	430	550	Loop	285 mm	400 mm	Slideout	10 kg	200 kg
mm	mm	mm	mm	mm	mm	Loop					
500	630	955	1005	430	550	Padded	215 mm	400 mm	Slidoout	11 kg	200 kg
mm	mm	mm	mm	mm	mm	Removable	וווווו כבכ	400 111111	Silueout	11 kg	200 kg
550	670	955	1005	460	550	Padded	315 mm	405 mm	Slideout	12 ka	300 kg
mm	mm	mm	mm	mm	mm	Removable	213 111111	ווווו כטדי	Silueout	12 Kg	Juu kg
600	720	955	1005	470	550	Padded	215 mm	405 mm	Slidoout	12 ka	300 kg
mm	mm	mm	mm	mm	mm	Removable	213 111111	ווווו כטדי	Silueout	13 kg	300 kg
650	770	955	1005	480		Padded	215 mm	405 mm	Slidoout	12 ka	300 kg
mm	mm	mm	mm	mm	mm	Removable	וווווו כדכ	ווווו כטד	Silueout	13 Kg	Juu ky