



KCARE TRANSIT MOBILE SHOWER COMMODE WITH SLIDEOUT FOOTPLATE



SKU: KA112

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Transit Mobile Shower Commode with Slideout Footplate is a versatile and practical solution for individuals with mobility difficulties. This commode features a slideout footplate that provides additional support for the user's feet and allows for easy side transfers. It is designed with a comfortable padded seat and a 3-strap backrest for optimal comfort and support. The commode is equipped with four 125mm total locking castors and an integrated push-handle for easy manoeuvrability. Additionally, it has swing-up arms that enable unrestricted side transfer and a removable pan for easy cleaning.

Please note: Padded seat is not included

FEATURES

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.
- 3 strap backrest for comfort and easy access.
- Integrated push-handle for manoeuvrability.
- Swing-up arms for unrestricted side transfer.
- Slide out/retractable non-weight bearing platform.
- 4x 125mm total locking castors.
- Designed to be used with a K Care patented removable seats.
- Available in 445, 500, 550, 600 and 650mm seat widths.
- **Seat, pan, and pan carrier not included - sold separately**

SPECIFICATIONS

Seat Width	Overall Width	Overall Depth	Overall Height	Seat Depth	Floor To Seat Height	Armrest Type	Armrest Height	Backrest Height	Footrest Type	Product Weight	Max User Weight
445 mm	480 mm	840 mm	1005 mm	430 mm	550 mm	Loop	285 mm	400 mm	Slideout	10 kg	200 kg
500 mm	630 mm	955 mm	1005 mm	430 mm	550 mm	Padded Removable	315 mm	400 mm	Slideout	11 kg	200 kg
550 mm	670 mm	955 mm	1005 mm	460 mm	550 mm	Padded Removable	315 mm	405 mm	Slideout	12 kg	300 kg
600 mm	720 mm	955 mm	1005 mm	470 mm	550 mm	Padded Removable	315 mm	405 mm	Slideout	13 kg	300 kg
650 mm	770 mm	955 mm	1005 mm	480 mm	550 mm	Padded Removable	315 mm	405 mm	Slideout	13 kg	300 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au