



## KCARE TRANSIT MOBILE SHOWER COMMODE WITH SWINGAWAY FOOTRESTS



**SKU:** KA114S-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The KCare Transit Mobile Shower Commode with Swing away Footrests is a versatile and durable product designed for individuals with mobility challenges. This commode can be used as a shower chair, toilet chair, or transport chair, providing comfort and convenience in one easy-to-use device.

The commode features swing away and removable footrests that can be adjusted to accommodate various leg lengths. The padded seat and backrest provide extra comfort, while the 3-strap backrest ensures secure positioning. The integrated push-handle allows for easy manoeuvrability, and the swing-up arms allow for unrestricted side transfer.

The commode is designed to be used with K Care patented removable seats and is available in various seat widths. It has four 125 mm total locking castors for added stability and safety.

### FEATURES

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.
- 3 strap backrest for comfort and easy access.
- Integrated push-handle for manoeuvrability.
- Swing-up arms for unrestricted side transfer.
- Swing away and removable height adjustable footrest.
- 4x 125mm total locking castors.
- Designed to be used with a K Care patented removable seats.
- **Seat, pan, and pan carrier not included - sold separately**

There is also a height adjustable version available, [click to view the KCare Height Adjustable Mobile](#)

[Shower Commode.](#)

## SPECIFICATIONS

Seat Width	Overall Depth	Overall Height	Seat Depth	Overall Width	Floor To Seat Height	Armrest Type	Armrest Height	Backrest Height	Footrest Type	Product Weight	Max User Weight
320 mm	750 mm	1005 mm	320 mm	485 mm	550 mm	Loop	285 mm	400 mm	Swing Away	10 kg	120 kg
400 mm	905 mm	1005 mm	430 mm	480 mm	550 mm	Loop	285 mm	400 mm	Swing Away	10 kg	200 kg
445 mm	890 mm	1005 mm	430 mm	520 mm	550 mm	Loop	285 mm	400 mm	Swing Away	10 kg	200 kg
500 mm	930 mm	1005 mm	430 mm	630 mm	550 mm	Padded Loop	315 mm	400 mm	Swing Away	11 kg	200 kg
550 mm	930 mm	1005 mm	460 mm	665 mm	550 mm	Padded Removable	315 mm	405 mm	Swing Away	12 kg	300 kg
600 mm	930 mm	1005 mm	470 mm	715 mm	550 mm	Padded Removable	315 mm	405 mm	Swing Away	13 kg	300 kg
650 mm	930 mm	1005 mm	480 mm	765 mm	550 mm	Padded Removable	315 mm	405 mm	Swing Away	14 kg	300 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)