



# KCARE TRANSIT MOBILE SHOWER COMMODE WITH SWINGAWAY FOOTRESTS

\$995.00 - \$2,126.00Price range: \$995.00 through \$2,126.00



**SKU:** KA114S-H

#### PRODUCT DESCRIPTION

#### **PURPOSE & BENEFITS**

The KCare Transit Mobile Shower Commode with Swing away Footrests is a versatile and durable product designed for individuals with mobility challenges. This commode can be used as a shower chair, toilet chair, or transport chair, providing comfort and convenience in one easy-to-use device.

The commode features swing away and removable footrests that can be adjusted to accommodate various leg lengths. The padded seat and backrest provide extra comfort, while the 3-strap backrest ensures secure positioning. The integrated push-handle allows for easy manoeuvrability, and the swing-up arms allow for unrestricted side transfer.

The commode is designed to be used with K Care patented removable seats and is available in various seat widths. It has four 125 mm total locking castors for added stability and safety.

#### **FEATURES**

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.
- 3 strap backrest for comfort and easy access.
- Integrated push-handle for manoeuvrability.
- Swing-up arms for unrestricted side transfer.
- Swing away and removable height adjustable footrest.
- 4x 125mm total locking castors.
- Designed to be used with a K Care patented removable seats.
- Seat, pan, and pan carrier not included sold separately

There is also a height adjustable version available, click to view the KCare Height Adjustable Mobile



## **Shower Commode.**

### **SPECIFICATIONS**

					Floor						Max
	: Overall h Depth				To Seat	Armrest Type	Armrest Height	Backrest Height	: Footrest Type	Product Weight	User
	•		•		Height						Weight
320	750	1005	320	485	550	Loop	285 mm	400 mm	Swing	10 kg	120 kg
mm	mm	mm	mm	mm	mm				Away		
400	905	1005	430	480	550	Loop	285 mm	400 mm	Swing	10 kg	200 kg
mm	mm	mm	mm	mm	mm				Away		
445	890	1005	430	520	550	Loop	285 mm	400 mm	Swing	10 kg	200 kg
mm	mm	mm	mm	mm	mm				Away	_ 0 9	_
500	930	1005	430	630	550	Padded 315	315 mm	400 mm	Swing	11 kg	200 kg
mm	mm	mm	mm	mm	mm	Loop	313		Away		
550	930	1005	460	665	550	Padded	315 mm	405 mm	Swing	12 kg	300 kg
mm	mm	mm	mm	mm	mm	Removable	)	.05 /////	Away	9	SOU NG
600	930	1005	470	715	550	Padded	315 mm	405 mm	Swing	13 kg	300 kg
mm	mm	mm	mm	mm	mm	Removable	313 111111	405 111111	Away	13 Kg	Joo kg
650	930	1005	480	765	550	Padded	315 mm	405 mm	Swing	14 kg	300 kg
mm	mm	mm	mm	mm	mm	Removable		100 111111	Away	± 1 Ng	200 119