



\*Seat not included

## KCARE TRANSIT MOBILE SHOWER COMMODE WITH WEIGHT BEARING FOOTPLATE

**\$1,105.00 - \$2,362.00** Price range: \$1,105.00 through \$2,362.00



**SKU: KA113**

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The KCare Transit Mobile Shower Commode with Weight Bearing Footplate is a versatile and practical solution for individuals with mobility limitations. This commode chair is designed for easy manoeuvrability, allowing it to be used in tight spaces or narrow hallways. The weight-bearing footplate provides added stability and safety for users, making it an ideal choice for those with reduced lower limb strength.

The chair is constructed with a sturdy, rust-resistant frame that can support up to 200-300 kg of weight. The seat and backrest are made from a comfortable, waterproof material that is easy to clean and maintain. The commode features a removable pan with lid for convenient and hygienic waste disposal.

This mobile shower commode chair is also equipped with slide-out footrests and loop-style armrests that can be easily adjusted to accommodate the user's specific needs. The chair's compact design makes it suitable for use in home care settings, nursing homes, hospitals, and other healthcare facilities.

**Please note: Padded seat not included**

### FEATURES

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.
- 3-strap backrest for comfort and easy access.
- Integrated push-handle for manoeuvrability.
- Swing-up arms for unrestricted side transfer.

- Slide out/retractable weight-bearing platform.
- 4x 125mm total locking castors.
- Designed to be used with a K Care patented removable seat.
- Available in 445mm and 500mm seat widths.
- **[Seat](#), [pan](#), and [pan carrier](#) not included - sold separately**

## SPECIFICATIONS

Seat Width	Overall Width	Overall Depth	Overall Height	Seat Depth	Floor To Seat Height	Armrest Type	Armrest Height	Backrest Height	Footrest Type	Product Weight	Max User Weight
445 mm	520 mm	830 mm	1005 mm	430 mm	550 mm	Loop	285 mm	400 mm	Weight Bearing	10 kg	200 kg
650 mm	766 mm	813 mm	1034 mm	480 mm	550 mm	Padded Removable	315 mm	405 mm	Weight Bearing	13 kg	300 kg
500 mm	630 mm	830 mm	1005 mm	430 mm	550 mm	Loop	285 mm	400 mm	Weight Bearing	11 kg	200 kg