



## KCARE TRAY MOBILE WALKER

**\$212.00**



**SKU:** KA354V

## PRODUCT DESCRIPTION

## PURPOSE & BENEFITS

This lightweight double tray mobile provides the user with greater independence. Ideal for in house use, the four castors allow for easy manoeuvrability on smooth surfaces.

## FEATURES

- Four castors for easy manoeuvrability
- Gives the user greater independence
- It comes with removable trays for easy cleaning
- Height can be adjusted by a push button system on each leg
- Removable trays for easy cleaning
- Ideal for home use and rehabilitation

## SPECIFICATIONS

Product Width (cm)	45
Product Depth (cm)	50
Product Height (cm)	79-92
Weight Capacity (kg)	100
DVA RAP Schedule	AH14
DVA RAP Code	KA354V
Handle Height Min (cm)	79

Handle Height Max (cm)

92