



KCARE TRAY MOBILE WALKER

\$212.00



SKU: KA354V

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

This lightweight double tray mobile provides the user with greater independence. Ideal for in house use, the four castors allow for easy manoeuvrability on smooth surfaces.

FEATURES

- Four castors for easy manoeuvrability
- Gives the user greater independence
- It comes with removable trays for easy cleaning
- Height can be adjusted by a push button system on each leg
- Removable trays for easy cleaning
- Ideal for home use and rehabilitation

SPECIFICATIONS

| | |
|------------------------|--------|
| Product Width (cm) | 45 |
| Product Depth (cm) | 50 |
| Product Height (cm) | 79-92 |
| Weight Capacity (kg) | 100 |
| DVA RAP Schedule | AH14 |
| DVA RAP Code | KA354V |
| Handle Height Min (cm) | 79 |

Handle Height Max (cm)

92