



KCARE ULTRA RYDER CARBON FIBRE FOLDING POWERCHAIR



SKU: KA390

PRODUCT DESCRIPTION

The **KCare Ultra Ryder Powerchair** is a compact, folding, carbon fibre power wheelchair designed for easy travel, storage, and daily use.

Features & Benefits:

- Versatile power wheelchair offers durability and ease of use
- Excellent manoeuvrability makes it perfect for indoor environments and local shops
- Smooth, effortless control for everyday activities
- Lightweight and compact front folding frame
- Swing back armrests for easier transfers
- Front suspension for a softer and smoother ride
- Remote control for attendant use
- Comfortable padded seat and back support
- Perfect for transport to family outings or holidays

Specifications:

- Overall Dimensions: 620mm W x 920mm D x 1000mm H
- Folded Dimensions: 620mm W x 770mm D x 320mm H
- Floor to Seat Height: 560mm
- Seat Size: 410mm W x 420mm D
- Width Between Armrests: 465mm
- Front Wheels: Aluminium, solid 180mm
- Rear Wheels: Magnesium, solid 275mm
- Net Weight: 17kg (excl. battery and accessories)
- Battery Weight: 1.8kg
- Battery Range: Minimum 20km*
- Top Speed: 6km/h

- Max User Weight: 150kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au