



KCARE ULTRA RYDER CARBON FIBRE FOLDING POWERCHAIR

\$4,025.00



SKU: KA390

PRODUCT DESCRIPTION

The **KCare Ultra Ryder Powerchair** is a compact, folding, carbon fibre power wheelchair designed for easy travel, storage, and daily use.

Features & Benefits:

- Versatile power wheelchair offers durability and ease of use
- Excellent manoeuvrability makes it perfect for indoor environments and local shops
- Smooth, effortless control for everyday activities
- Lightweight and compact front folding frame
- Swing back armrests for easier transfers
- Front suspension for a softer and smoother ride
- · Remote control for attendant use
- Comfortable padded seat and back support
- · Perfect for transport to family outings or holidays

Specifications:

Overall Dimensions: 620mm W x 920mm D x 1000mm H
 Folded Dimensions: 620mm W x 770mm D x 320mm H

• Floor to Seat Height: 560mm

Seat Size: 410mm W x 420mm D
Width Between Armrests: 465mm

Front Wheels: Aluminium, solid 180mm
Rear Wheels: Magnesium, solid 275mm

Net Weight: 17kg (excl. battery and accessories)

Battery Weight: 1.8kg

• Battery Range: Minimum 20km*



• Top Speed: 6km/h

• Max User Weight: 150kg