



LIGHTWEIGHT KITCHEN PERCHING STOOL - HEIGHT ADJUSTABLE



SKU: HOM-549B

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

This perching stool is designed to take the weight off your feet and back to ease pain and reduce fatigue when you are working at a benchtop for example in preparing meals. The stool has a broad sloping padded white vinyl seat for maximum comfort and has large non-slip feet to keep it stable when you have your weight leaning on it.

It is designed to be sat on with the angle of the seat ensuring that weight is shared between your feet and the seat so you are "perching" and can easily move to and away from the seat for pause breaks whilst involved in preparation tasks.

Adjustable height with arms, padded seat and padded backrest.

SPECIFICATIONS

- Width between arms: 455mm
- Seat depth: 300mm
- Seat to top of backrest 300mm
- Leg height adjustment floor to seat: 500mm - 655mm
- User weight up to 100kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au