



OPPO KNEE IMMOBILISER LENGTH 18"

~~\$57.00~~ Original price was: ~~\$57.00~~. ~~\$42.00~~ Current price is: \$42.00. While Stocks Last Online Only



SKU: OPP4030-1860

PRODUCT DESCRIPTION

Purpose & Benefits

The **OPPO Knee Immobiliser** (18-inch length, XL) provides strong yet comfortable immobilisation for the knee joint, particularly suited for post-injury or post-operative recovery. Its open-front design and adjustable upper and lower straps allow for simple application and a tailored fit, even for limited mobility users. Designed with multiple aluminium stays to maintain alignment, it minimises unwanted movement and includes a pressure-relief opening around the kneecap to enhance comfort. This support encourages a stable environment conducive to healing, offering reassurance and reducing the risk of further injury during recovery.

Features

- Aluminium stays along each side provide rigid structural support and restrict knee movement
- Upper and lower adjustable straps deliver a customised, secure fit
- Open-front design for effortless application and adaptability
- Patella relief zone enhances comfort by reducing pressure on the kneecap

Specifications

- **Brand:** OPPO Medical

- **Size:** XL
- **Length:** 18 in (45.7 cm)
- **Material Composition:** Approx. 80 % polyester, 15 % foam, 5 % nylon