



KNIGHT COMFORT BARIATRIC CHAIR

\$661.00



SKU: AGE400B

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Australian-made Knight Comfort Bariatric Chair is a premium quality design engineered for comfort. It is a fully upholstered, multi-purpose, support chair that is perfect for use in any environment where an adjustable floor-to-seat height is required, such as the lounge, dining room, or bedroom.

The Knight Comfort Chair is available in a variety of colour options in both vinyl and fabric.

FEATURES

- Knight Bariatric Healthcare Chair
- Wider Seat
- Custom Designed Heavy-Duty Metal Frame
- Frames are Finished in a Durable Powder-coat
- Seat & Back Rests are Solid Ply Timber
- Foams are Dunlop Enduro Fire Retardant
- Padded Armrests
- Excellent Comfort & Support
- Commercial Grade Upholstery
- High-Quality Stitching
- Hand Built in Australia
- Floor to Seat Adjustment
- 350kg Weight Capacity
- 10-Year Warranty

COLOUR & STYLE

Some colours in the image have been discontinued, select a colour from the available list of options.

DIMENSIONS

Bariatric

Width	650mm
Depth	450mm
Floor to Seat Adjustment	460mm-585mm
Back Rest Height	490mm
Arm Rest Height	180mm
Weight Capacity	400kg
Package Weight	15kg

Buy the Knight Comfort Bariatric Chair with the [NDIS payment method](#)