




KNIGHT COMFORT BARIATRIC CHAIR



SKU: AGE400B

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Australian-made Knight Comfort Bariatric Chair is a premium quality design engineered for comfort. It is a fully upholstered, multi-purpose, support chair that is perfect for use in any environment where an adjustable floor-to-seat height is required, such as the lounge, dining room, or bedroom.

The Knight Comfort Chair is available in a variety of colour options in both vinyl and fabric.

FEATURES

- Knight Bariatric Healthcare Chair
- Wider Seat
- Custom Designed Heavy-Duty Metal Frame
- Frames are Finished in a Durable Powder-coat
- Seat & Back Rests are Solid Ply Timber
- Foams are Dunlop Enduro Fire Retardant
- Padded Armrests
- Excellent Comfort & Support
- Commercial Grade Upholstery
- High-Quality Stitching
- Hand Built in Australia
- Floor to Seat Adjustment
- 350kg Weight Capacity
- 10-Year Warranty

COLOUR & STYLE

Some colours in the image have been discontinued, select a colour from the available list of options.

DIMENSIONS

Bariatric

Width	650mm
Depth	450mm
Floor to Seat Adjustment	460mm-585mm
Back Rest Height	490mm
Arm Rest Height	180mm
Weight Capacity	400kg
Package Weight	15kg

Buy the Knight Comfort Bariatric Chair with the [NDIS payment method](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au