



KNIGHT COMFORT STANDARD CHAIR

\$506.00 Original price was: \$506.00.<u>\$455.00</u>Current price is: \$455.00.



SKU: AGE400

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Australian-made Knight Comfort chair is a premium quality design engineered for comfort.

A fully upholstered, multi-purpose, support chair is designed for use in any environment where an adjustable floor to seat height is required: Lounge, dining, bedroom.

Available in a variety of colour options in both vinyl and fabric.

FEATURES

- Knight Standard Healthcare Chair
- Wider Seat
- Custom Designed Heavy-Duty Metal Frame
- Frames are Finished in a Durable Powder-coat
- Seat & Back Rests are Solid Ply Timber
- Foams are Dunlop Enduro Fire Retardant
- Padded Armrests
- Excellent Comfort & Support
- Commercial Grade Upholstery
- High-Quality Stitching
- Hand Built in Australia
- Floor to Seat Adjustment
- 350kg Weight Capacity
- 10-Year Warranty



COLOUR & STYLE

Some colours in the image have been discontinued, select a colour from the available list of options.

DIMENSIONS

Standard

Width	450mm
Depth	450mm
Floor to Seat Adjustment	460mm-585mm
Back Rest Height	490mm
Arm Rest Height	180mm
Weight Capacity	350kg
Package Weight	14kg