



## **KNIGHT COMFORT STANDARD CHAIR**

**\$506.00**



**SKU: AGE400**

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The Australian-made Knight Comfort chair is a premium quality design engineered for comfort.

A fully upholstered, multi-purpose, support chair is designed for use in any environment where an adjustable floor to seat height is required: Lounge, dining, bedroom.

Available in a variety of colour options in both vinyl and fabric.

### **FEATURES**

- Knight Standard Healthcare Chair
- Wider Seat
- Custom Designed Heavy-Duty Metal Frame
- Frames are Finished in a Durable Powder-coat
- Seat & Back Rests are Solid Ply Timber
- Foams are Dunlop Enduro Fire Retardant
- Padded Armrests
- Excellent Comfort & Support
- Commercial Grade Upholstery
- High-Quality Stitching
- Hand Built in Australia
- Floor to Seat Adjustment
- 350kg Weight Capacity
- 10-Year Warranty

## COLOUR & STYLE

Some colours in the image have been discontinued, select a colour from the available list of options.

## DIMENSIONS

### Standard

Width	450mm
Depth	450mm
Floor to Seat Adjustment	460mm-585mm
Back Rest Height	490mm
Arm Rest Height	180mm
Weight Capacity	350kg
Package Weight	14kg