

LANFRANCO BONDI RECLINER



SKU: LANBON

PRODUCT DESCRIPTION

Options

The Bondi can be supplied in a small, medium, or large size, and with either a single or dual motor variation.

All Lanfranco products are available in a range of colours and fabrics. Many can be found by clicking the link below.

[Click here to view the available colour range.](#)

Specifications

Frame Type: Frames incorporate Radiata Plantation Timbers, Australia's renewable resource.

Suspension Warranty: 10 year warranty - Black Cat Webbing.

Frame Warranty: Lifetime warranty.

Recliner Release: Handle or trigger action.

Back Cushion: Top grade silicon coated fibre - removable.

Seat Cushion: Dunlop Enduro - lifetime warranty.

FINISHED OVERALL DIMENSIONS:

	SMALL LIFT RECLINER	MEDIUM LIFT CHAIR	LARGE LIFT CHAIR
HEIGHT	overall 1060mm / 530mm floor to seat	overall 1100mm / 530mm floor to seat	overall 1140mm / 530mm floor to seat
WIDTH	overall 780mm / 400mm between arms	overall 810mm / 450mm between arms	overall 870mm / 510mm between arms
DEPTH	overall 930mm / 500mm seat depth	overall 930mm / 500mm seat depth	overall 930mm / 500mm seat depth

WARRANTY

2 Years Warranty: Mechanism, fabric and electronics

10 Years Warranty: Timber Frame and Foam.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au

