




LEVABO TURN ALL BARIATRIC BED



SKU: PCLV-70098

PRODUCT DESCRIPTION

Purpose & Benefits

The **Levabo Turn All Bariatric Bed** is an advanced pressure care and positioning system designed to support the needs of bariatric patients. It provides automated lateral turning to reduce the risk of pressure injuries, assist with pulmonary care, and improve circulation. Designed for comfort and efficiency, the Turn All system enhances patient safety while reducing physical strain on caregivers in clinical and home care environments.

Features

- Automated lateral turning system to assist with pressure relief and repositioning
- Designed for bariatric users, providing stable and even weight distribution
- Supports prevention and management of pressure ulcers by promoting regular movement
- Programmable turn cycles to adjust the degree and timing of rotation
- Reduces caregiver workload by automating patient repositioning routines
- Comfortable, supportive surface designed for long-term clinical use
- Easy integration into existing hospital or home care bed setups

Specifications

- **Model:** Levabo Turn All Bariatric Bed
- **Function:** Automated lateral turning for pressure care and repositioning
- **User Type:** Bariatric patients requiring pressure injury prevention and mobility assistance
- **Turn Cycle:** Adjustable rotation intervals and angles
- **Care Setting:** Suitable for hospital, aged care and home environments