



LIFTSEAT HOME PLUS



SKU: LSH2

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

LiftSeat Home Plus is a powered toilet lift that helps people with mobility challenges get on and off the toilet safely and easily. It is a free-standing unit that fits over any standard toilet, and it is powered by a rechargeable battery. The LiftSeat Home Plus has a weight capacity of up to 203 kg, making it suitable for bariatric users.

The LiftSeat Home Plus is a great way to improve your independence and reduce your risk of falls. It is also a valuable tool for caregivers, as it can help them to provide safer and more comfortable toileting assistance.

FEATURES

- Powered toilet lift
- Free-standing unit
- Fits over any standard toilet
- Rechargeable battery
- Weight capacity of up to 203 kg
- Built-in safety features
- Easy to use
- Ideal for people with mobility challenges and caregivers

If you are looking for a safe and easy way to get on and off the toilet, the LiftSeat Home Plus is a great option for you. It is a reliable and durable product that can help you to maintain your independence and dignity.

SPECIFICATIONS

- Weight capacity: 203 kg
- Lift height: 500 mm
- Seat width: 460 mm
- Seat depth: 400 mm
- Seat height: 460 mm
- Base dimensions: 700 x 520 mm
- Battery life: 120 lifts
- Charging time: 4 hours

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au