



LIFTSEAT HOME



SKU: LSH1

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

LiftSeat Home is a powered toilet lift that helps people with mobility challenges toilet independently or with reduced assistance. It is a safe and easy-to-use device that lowers the user onto the toilet and raises them back up again. The LiftSeat Home is powered by a quiet motor and has a smooth, controlled movement that minimizes the risk of falls.

The LiftSeat Home is easy to install and can be used over any standard toilet. It is also adjustable to fit most users, and it has a weight capacity of 146 kg.

FEATURES

- Powered by a quiet motor for smooth, controlled movement
- Adjustable to fit most users
- Weight capacity of 146 kg
- Easy to install
- Can be used over any standard toilet
- Reduces the risk of falls and injuries
- Promotes independence and dignity

The LiftSeat Home is a valuable tool for people with mobility challenges. It can help them regain their independence and dignity, and it can also help to reduce the risk of falls and injuries.

SPECIFICATIONS

- Weight capacity: 146 kg
- Dimensions: 68.58 x 54.61 x 43.18 cm

- Weight: 24kg
- Power: 120V AC, 50/60Hz
- Battery backup: 12V DC, 7 amp-hour
- Operation time: Up to 20 lifts on battery power

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au