



LIFTSEAT VERTICA

\$8,250.00



SKU: LSV1

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

LiftSeat Vertica is a powered toilet lift that helps people with mobility challenges toilet independently and safely. It rises vertically upwards with no tilt at any point in its movement, making it much easier to lock the knees before standing. The LiftSeat Vertica is also easy to use, with a simple 2-button handheld control.

FEATURES

- Vertical lift with no tilt
- Smooth and easy operation
- Simple 2-button handheld control
- Compact and lightweight design
- Safe working load of 203kg

SPECIFICATIONS

- Overall dimensions: 700 x 450 x 1000 mm
- Seat dimensions: 400 x 400 mm
- Weight: 25 kg
- Safe working load: 200 kg
- Lifting range: 400 mm
- Power supply: 230 V AC
- Power consumption: 150 W