



LIFTSEAT VERTICA



SKU: LSV1

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

LiftSeat Vertica is a powered toilet lift that helps people with mobility challenges toilet independently and safely. It rises vertically upwards with no tilt at any point in its movement, making it much easier to lock the knees before standing. The LiftSeat Vertica is also easy to use, with a simple 2-button handheld control.

FEATURES

- Vertical lift with no tilt
- Smooth and easy operation
- Simple 2-button handheld control
- · Compact and lightweight design
- Safe working load of 203kg

SPECIFICATIONS

• Overall dimensions: 700 x 450 x 1000 mm

• Seat dimensions: 400 x 400 mm

• Weight: 25 kg

Safe working load: 200 kg
Lifting range: 400 mm
Power supply: 230 V AC
Power consumption: 150 W



For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773 brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au