



## **LIGHT DRIVE RIGID WHEELCHAIR POWER ASSIST | TURN ANY WHEELCHAIR INTO A POWER CHAIR**



**SKU:** 81M21CF

### **PRODUCT DESCRIPTION**

#### **PURPOSE & BENEFITS**

The Light Drive power assist is an electric motor that turns your manual wheelchair into an ultra-light, high-performance electric wheelchair in seconds!

An extended colour range, a new anti-tip, and a more ergonomic clutch system. The Light Drive power assist is changing in order to make life even more beautiful.

The motor attachment simply hooks onto the rear of your rigid or folding manual wheelchair.

The Light Drive unit is operated by a joystick while maintaining the chair's manoeuvrability – it stops when you let go of the joystick and has a 360-degree movement including full reversal.

By releasing a clutch, you can free the wheel giving you the choice between using the power add-on assistance or manual use.

Safety is an important factor of the Light Drive with a downhill braking feature in place that automatically controls your descent.

Benoit Systemes have designed the Light Drive Wheelchair Power Assist to be a robust & reliable power add for most manual wheelchairs.

#### **FEATURES**

- Turns your self-propelled wheelchair into a power chair in seconds
- Climb hills with no effort, coast across the grass, and sweep along gravel paths and uneven pavers!

- Top speed 10km/h
- Range 14km-16km on one battery charge
- Funded by NDIS
- The unit is easy to unhook from the wheelchair making it an in-expensive power device to transport via a vehicle – no need for an expensive adaptive van for transportation.
- Provides wheelers with extra power while maintaining exercise, muscle tone, strength and energy.
- Avoids the need to switch to a heavy power chair.
- Can be fitted to kids' wheelchairs with a minimum width of 40 cm between wheels.

## SPECIFICATIONS

Max. Weight	Approx. 100 kg (16 stone) for optimum performance. Can go up to 120 kg (9 stone) - terrain dependent.
Weight of Unit	Motor 9kg. Battery 5kg. Brackets & Joystick 500g.
Climb	20% (12°) for a person weighing 70kg (11 stone)
Speed	Adjustable between 0-10 km/h (0-6 mph)
Range	14 to 16 km
Tyres	Contact us for recommendations. You need tyres with excellent traction and a very thick layer of puncture-resistant material. Solid tyres are not compatible.
Front Castors	No smaller than 6" in diameter.
Battery	Lithium
Control	Joystick (programmable)

<https://youtu.be/kPQmB3Tyg-E>

Buy this Light Drive Wheelchair Power Assist online at Medimart using Stripe, PayPal, Afterpay or [NDIS method](#) today!

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

**SA**

**Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

**Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

**Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

**VIC**

**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

**Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

**NT**

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)