



## LOCOMO SALJOL ALUMINIUM WALKER

**\$1,221.00**



**SKU: 42101-H**

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Introducing the Locomo Saljol Aluminium Walker, an exceptional mobility solution designed with a perfect blend of strength, versatility, and convenience. Constructed from high-quality aluminium, this walker ensures durability while maintaining a lightweight profile, making it ideal for individuals seeking support and stability during their daily activities.

### FEATURES

- **Aluminium Frame:** The walker is crafted from premium aluminium, guaranteeing robustness while ensuring it remains light enough for easy maneuverability.
- **Adjustable Height:** Users can tailor the height of the walker to suit their individual needs, providing a comfortable and ergonomic experience.
- **Ergonomic Handles:** The handles are designed to offer maximum comfort, reducing strain on the wrists and hands and allowing for prolonged use without discomfort.
- **Compact Folding Design:** Storage and transportation are made easy with the walker's compact folding design, ensuring it can be stowed away or transported with minimal hassle.
- **Sleek Aesthetic:** Sporting a modern and stylish look, this walker combines functionality with elegance, enhancing the user's mobility aid collection.
- **Durable Wheels:** Designed to tackle various terrains, the durable wheels offer a smooth and stable ride, enhancing mobility and confidence.
- **Reliable Brake System:** Safety is a top priority, and the walker features a reliable brake system, allowing users to securely lock it in place when needed.

## **SPECIFICATIONS**

- Rollator weight: 6.6kg
- Maximum user weight capacity: 150kg
- Maximum load of the bag: 5kg
- Seat height: 62cm
- Handle height range: 81 – 103cm
- Seat width: 46cm
- Overall width: 61cm
- Width folded: 21cm
- Front tyres: 250mm