




LOCOMO SALJOL CARBON WALKER

\$1,395.00 - \$1,495.00 Price range: \$1,395.00 through \$1,495.00



SKU: 41103-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Experience unparalleled mobility and style with the Locomo Saljol Carbon Walker, a masterclass in innovation and ergonomic design. Crafted from high-quality carbon fibre, this walker sets a new standard in lightweight, durable mobility aids, ensuring users can navigate their environment with ease and confidence.

FEATURES

- **Carbon Fiber Construction:** The Locomo Saljol is built with carbon fibre, renowned for its strength-to-weight ratio, ensuring the walker is not just lightweight but also exceptionally sturdy.
- **Adjustable Height:** Catering to individual user needs, the walker's height can be easily adjusted, providing a tailored fit for optimal comfort and support.
- **Ergonomic Handles:** Designed with user comfort in mind, the ergonomic handles reduce strain on the wrists and hands, allowing for prolonged use without discomfort.
- **Foldable Design:** Convenience is key, and the Locomo Saljol delivers with a foldable design, making storage and transportation hassle-free.
- **Stylish Aesthetic:** With a sleek and modern look, this walker not only provides functionality but also adds a touch of elegance to your mobility aid collection.
- **Brake System:** Safety is paramount, and the integrated brake system ensures users can secure the walker in place as needed.
- **Accessory Options:** Customise your walker with the included bag and cane holder

Discover the perfect balance of style, strength, and functionality with the Locomo Saljol Carbon Walker. Navigate your world with confidence, knowing you have a top-tier mobility aid at your

service. Elevate your mobility experience and step into a world of unmatched elegance and ergonomic design with the Locomo Saljol.

SPECIFICATIONS

- Rollator weight: 5.6kg (*without the bag, 6kgs with the bag*)
- User height: >160cm
- Maximum weight capacity: 150kg
- Maximum load of the bag: 5kg
- Seat height: 62cm
- Handle height range: 81 – 103cm
- Seat width: 45cm
- Overall width: 61cm
- Width folded: 21cm

ACCESSORIES

- [Saljol Walker Backrest](#)