




## LUXA HOMECARE BED



**SKU:** LUXA-P

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Luxa represents the apex of Homecare bed design, offering exceptional comfort and care, all the while preserving an individual's dignity and augmenting the aesthetics of their bedroom.

This bed comes equipped with a multitude of high-quality features. These include a backrest with auto-regression capabilities, a manual vascular support system, and an adjustable height ranging from a lofty 725mm to a mere 100mm. Additionally, it provides the possibility to integrate Bluetooth handsets.

Moreover, it is capable of being upgraded with several optional accessories. Among these are a robust horizontal rail, provisions for the inclusion of self-help poles, lights underneath the bed, and many more.

### FEATURES

- Extremely Low Minimum Position
- Vertical Lift
- Auto-Regression Backrest
- Ash or Ocean Padded Headboard with full
- Padded Surround

### SPECIFICATIONS

Bed Type	Overall Dimensions (mm)	Sleeping Surface (mm)	SWL (kg)	Low/High Position (mm)	Back Angle (degrees)	Knee Angle (degrees)	Reverse Angle (degrees)	Trendelenburg Angle (degrees)
----------	-------------------------	-----------------------	----------	------------------------	----------------------	----------------------	-------------------------	-------------------------------

Bed Type	Overall Dimensions (mm)	Sleeping Surface (mm)	SWL (kg)	Low/High Position (mm)	Back Angle (degrees)	Knee Angle (degrees)	Reverse Angle (degrees)	Trendelenburg Angle (degrees)
Luxa Homecare	2160 x 1000	2000 x 900	250	150/725	62	35	15	15
Luxa King Single	2160 x 1170	2000 x 1070	250	150/725	62	35	15	15
Luxa Double	2160 x 1470	2000 x 1370	250	170/725	62	35	15	15
Luxa Queen	2160 x 1600	2000 x 1500	250	170/725	62	35	15	15

## BROCHURE

Click here to view the brochure [Luxa-by-Caremed-Alrick](#)