



## LUXA HOMECARE TIMBER FRAME BED



SKU: N/A

## PRODUCT DESCRIPTION

### Purpose & Benefits of Luxa Homecare Timber Frame Bed

The Luxa Homecare Timber Frame Bed is designed to provide a comfortable and safe sleeping environment for individuals in need of home care. This bed offers several features that promote convenience, accessibility, and durability, making it an ideal choice for those requiring long-term or temporary bed rest.

### Features

- **Adjustable Height:** The bed's height can be easily adjusted to accommodate different user needs, making it easier for caregivers to assist with tasks such as bed transfers and personal care.
- **Timber Frame:** Constructed from high-quality timber, the frame is sturdy and durable, providing a stable sleeping surface for users.
- **Head and Foot Adjustment:** Allows for customised positioning to enhance comfort while reading, watching TV, or sleeping.
- **Side Rails:** Provides added safety by preventing falls and offering support for users getting in and out of bed.
- **Easy Assembly:** Designed for simple and quick assembly without the need for special tools.
- **Mattress Options:** Compatible with a variety of mattresses to suit individual comfort preferences and medical needs.

## SPECIFICATIONS

Bed Type	Overall Dimensions (mm)	Sleeping Surface (mm)	SWL (kg)	Low/High Position (mm)	Back Angle (degrees)	Knee Angle (degrees)	Reverse Angle (degrees)	Trendelenburg Angle (degrees)
Luxa Homecare	2160 x 1000	2000 x 900	250	150/725	62	35	15	15
Luxa King Single	2160 x 1170	2000 x 1070	250	150/725	62	35	15	15
Luxa Double	2160 x 1470	2000 x 1370	250	170/725	62	35	15	15
Luxa Queen	2160 x 1600	2000 x 1500	250	170/725	62	35	15	15

## BROCHURE

Click here to view the brochure [Luxa-by-Caremed-Alrick](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriiffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)