




MEDICAL INFRARED FOREHEAD THERMOMETER



SKU: XL-F02

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Using medical infrared quality technology, this non-contact device ensures quick & easy results for both human & objects. The LCD screen displays readings clearly, with a built-in memory function for effortless tracking and multiple coloured displays to easily indicate if temperature are inside normal ranges.

- Easily & safely measure the whole family in seconds, with no need to clean between measurements
- Keep your workplace safe with quick effortless screening of all staff and visitors
- Check the temperature of food, and other objects to keep you and your environment safe

FEATURES

- Infrared Digital Thermometer
- Contactless Thermometer
- Fast, Accurate, Dependable Readings
- Easy to Read LED Screen
- 2 x AAA Batteries Included
- Auto Shut Off System in 20 seconds
- Memory Recall for 32 previous readings
- Convenient and Essential Item
- Ideal for use in Offices, Home, Hospitals, Doctors Surgery
- Travel Handy

SPECIFICATIONS

Product name:	Medical infrared thermometer
---------------	------------------------------

Model:	XL-F02
Display accuracy:	0.1°C / °F
Measurement range:	Person mode: 32.0°C ~ 42.9°C
Recording:	32 sets
Unit:	°F / °C
Period of use:	3 years
Measure time:	1 second
Function distance:	5 cm ~ 10 cm
Over time:	20 seconds
Relative humidity:	Less than or equal to 80%
Battery:	DC 3V (2 x AAA batteries)
Maximum error:	Within 35.0°C ~ 42.0°C $\pm 0.2^{\circ}\text{C}$; otherwise $\pm 0.3^{\circ}\text{C}$
Measurement (mm):	155mm x 155mm x 45mm (L x W x H)
Weight:	90g

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au