



## MEDICAL INFRARED FOREHEAD THERMOMETER

**\$99.95**



**SKU: XL-F02**

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Using medical infrared quality technology, this non-contact device ensures quick & easy results for both human & objects. The LCD screen displays readings clearly, with a built-in memory function for effortless tracking and multiple coloured displays to easily indicate if temperature are inside normal ranges.

- Easily & safely measure the whole family in seconds, with no need to clean between measurements
- Keep your workplace safe with quick effortless screening of all staff and visitors
- Check the temperature of food, and other objects to keep you and your environment safe

### FEATURES

- Infrared Digital Thermometer
- Contactless Thermometer
- Fast, Accurate, Dependable Readings
- Easy to Read LED Screen
- 2 x AAA Batteries Included
- Auto Shut Off System in 20 seconds
- Memory Recall for 32 previous readings
- Convenient and Essential Item
- Ideal for use in Offices, Home, Hospitals, Doctors Surgery
- Travel Handy

### SPECIFICATIONS

Product name:	Medical infrared thermometer
---------------	------------------------------

Model:	XL-F02
Display accuracy:	0.1°C / °F
Measurement range:	Person mode: 32.0°C ~ 42.9°C
Recording:	32 sets
Unit:	°F / °C
Period of use:	3 years
Measure time:	1 second
Function distance:	5 cm ~ 10 cm
Over time:	20 seconds
Relative humidity:	Less than or equal to 80%
Battery:	DC 3V (2 x AAA batteries)
Maximum error:	Within 35.0°C ~ 42.0°C $\pm 0.2^{\circ}\text{C}$ ; otherwise $\pm 0.3^{\circ}\text{C}$
Measurement (mm):	155mm x 155mm x 45mm (L x W x H)
Weight:	90g