

MEDIFONE SAFETY PENDANT



SKU: MEDIFONE

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The MediFone Safety Pendant is a state-of-the-art personal emergency response system designed for individuals who need an easy and reliable way to call for help in case of an emergency. This device is particularly beneficial for seniors, people with health issues, or those living alone, providing peace of mind for both users and their loved ones.

FEATURES

- **One-Button SOS:** Features a simple, large, easy-to-press SOS button that, when activated, sends an alert to pre-set emergency contacts or a monitoring centre.
- **Two-Way Communication:** Equipped with a built-in speaker and microphone, allowing for two-way voice communication with emergency responders or family members directly through the pendant.
- **GPS Location Tracking:** Includes GPS technology to pinpoint the user's location, ensuring that help can be sent to the exact location in a timely manner.
- **Fall Detection Sensor:** Some models may include automatic fall detection, which triggers an alert if a fall is detected, even if the user is unable to press the button.
- **Water-Resistant Design:** Designed to be water-resistant, making it safe to use in the shower or in humid conditions.
- **Long Battery Life:** Offers a long battery life with easy charging options, minimizing the need for frequent recharging.

- **Compact and Wearable:** Lightweight and can be worn around the neck, providing convenience and ensuring it is always within reach.

The MediFone Safety Pendant is an essential device for anyone who may be at risk of medical emergencies or falls, offering immediate access to assistance at the press of a button. It's a practical solution for maintaining independence while ensuring safety and security in various living environments.

SPECIFICATIONS

Size in millimetres / Weight	65 X 48 X 19 / 56 grams
The SOS button will contact	Up to 8 numbers – Text and/or Call
Sequential SOS calls go to each Emergency Contact	YES
Programmable to call landline, mobile and 000	YES
GPS location is provided by pressing the SOS or by:	Call and/or Text “123456F”
Programmable by:	SMS Text
Phone Service Carriers	Telstra’s 3G & 4G frequency bands
Available Speed Dial N°s	8
Capable of Receiving Calls	YES
Carrying Methods	Lanyard, Belt Clip
Recharges in a docking cradle	YES
Fall Detector	Yes. Cancellable within 15 seconds
Medicine Reminder	YES
Optional Accessories	Button Belt Pouch, Velcro Belt Pouch, Shower Pouch

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au