


MERITS 745 ECO SCOOTER



SKU: S745-1YL-H



PRODUCT DESCRIPTION

Purpose & Benefits of Merits 745 Eco Scooter

The **Merits 745 Eco Scooter** is a compact and practical mobility scooter designed for users who need a reliable and easy-to-use solution for everyday travel. With its user-friendly controls, streamlined design, and comfortable seating, the 745 Eco is ideal for short trips, errands, and indoor or light outdoor use.

This scooter offers excellent manoeuvrability in tight spaces, making it well-suited for shopping centres, walkways, and home environments. Its lightweight build and simple operation provide confidence and independence for users seeking an affordable and efficient mobility option.

Key Features:

- Compact 4-wheel design for stability and manoeuvrability
- Comfortable padded seat with flip-up armrests
- Simple control panel for intuitive operation
- Lightweight and easy to transport

- Suitable for indoor and light outdoor use
- Front basket for everyday convenience
- Ideal for users seeking reliable mobility on a budget

Specifications:

- Overall Dimension: 120 x 61 x 99cm
- Seat Dimension: 46 x 40 x 46cm
- Seat Height from Ground: 53-58cm
- Seat Height from Deck: 40-47cm
- Wheel Base: 85cm
- Max Speed 10kph
- Range 20-30km
- Turn Radius 157cm
- Ground clearance 9.5cm
- Motor DC 24V, 210W (Rating)/ 715W (Max)
- Controller S-Drive 70A
- Battery 12V/34A/H x 2pcs
- Charger 5A off-board Gradient 6°
- Front Wheel 10" Pneumatic tire
- Rear Wheel 10" Pneumatic tire

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au