




### **MERITS 745 ECO SCOOTER**

**\$3,299.24 - \$3,339.05** Price range: \$3,299.24 through \$3,339.05



**SKU:** S745-1YL-H

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits of Merits 745 Eco Scooter**

The **Merits 745 Eco Scooter** is a compact and practical mobility scooter designed for users who need a reliable and easy-to-use solution for everyday travel. With its user-friendly controls, streamlined design, and comfortable seating, the 745 Eco is ideal for short trips, errands, and indoor or light outdoor use.

This scooter offers excellent manoeuvrability in tight spaces, making it well-suited for shopping centres, walkways, and home environments. Its lightweight build and simple operation provide confidence and independence for users seeking an affordable and efficient mobility option.

### **Key Features:**

- Compact 4-wheel design for stability and manoeuvrability
- Comfortable padded seat with flip-up armrests
- Simple control panel for intuitive operation
- Lightweight and easy to transport

- Suitable for indoor and light outdoor use
- Front basket for everyday convenience
- Ideal for users seeking reliable mobility on a budget

## **Specifications:**

- Overall Dimension: 120 x 61 x 99cm
- Seat Dimension: 46 x 40 x 46cm
- Seat Height from Ground: 53-58cm
- Seat Height from Deck: 40-47cm
- Wheel Base: 85cm
- Max Speed 10kph
- Range 20-30km
- Turn Radius 157cm
- Ground clearance 9.5cm
- Motor DC 24V, 210W (Rating)/ 715W (Max)
- Controller S-Drive 70A
- Battery 12V/34A/H x 2pcs
- Charger 5A off-board Gradient 6°
- Front Wheel 10" Pneumatic tire
- Rear Wheel 10" Pneumatic tire