



MERITS 745 PLUS SCOOTER

\$4,430.60



SKU: S745A-10B-H

PRODUCT DESCRIPTION

Purpose & Benefits of Merits 745 Plus Scooter

The **Merits 745 Plus Scooter** is an upgraded compact mobility scooter designed to provide enhanced comfort, stability, and ease of use for individuals seeking a reliable solution for everyday travel. Built with a robust 4-wheel design and improved features over the standard model, the 745 Plus is ideal for both indoor environments and smooth outdoor surfaces.

It offers a comfortable padded seat, simple controls, and an easy disassembly system for convenient transport and storage. The 745 Plus is perfect for users looking for a portable and efficient scooter that doesn't compromise on quality or safety.

Key Features

- 4-wheel configuration for improved stability and control
- Comfortable swivel seat with padded armrests
- Easy disassembly for transport and storage
- Simple, intuitive tiller controls
- Suitable for shopping, appointments, and local outings

- Front basket for daily convenience
- Ideal for users who require a step above basic mobility scooters

Specifications

Total Size A x B x C	1320mm x 610mm x 1020mm
Seat Size D x E x F	460mm x 400mm x 460mm
Seat Height from Ground G	610mm - 660mm
Seat Height from Deck H	440mm - 490mm
Wheel Base I	980mm
Maximum Speed	10kph
Range	30km - 40km
Weight with Battery	98kg
Battery Weight	16kg
SWL	150kg