




MERITS AURORA 2 HILL CLIMBER 3 WHEEL SCOOTER



SKU: S836C-HC-BR-H

PRODUCT DESCRIPTION

Purpose & Benefits of Merits Aurora 2 Hill Climber 3 Wheel Scooter

The **Merits Aurora 2 Hill Climber Scooter** is a powerful and compact 3-wheel mobility scooter designed for users who require excellent hill-climbing capability and superior manoeuvrability. Ideal for both indoor use and outdoor travel on sloped terrain, the Aurora 2 provides smooth, stable performance with its robust motor and smart design.

Its 3-wheel configuration allows for a tighter turning radius, making it easy to navigate in shopping centres and narrow pathways. A padded high-back seat, intuitive controls, and front suspension ensure a comfortable ride, while the easy disassembly system makes it a practical choice for everyday mobility.

Key Features

- High-torque motor for enhanced hill-climbing performance
- 3-wheel design for tight turning and easy navigation in confined spaces
- Padded swivel seat with adjustable armrests for comfort and access
- Front suspension for a smoother ride over uneven surfaces

- Disassembles easily for transport and storage
- Ideal for shopping, local errands, and hilly outdoor environments
- Simple, user-friendly tiller controls
- Front basket for carrying personal items

Specifications

- Weight capacity: 182 kg
- Max speed: 12 kph (limited in Au to 10 km/h)
- Max speed (Hill Climber): 10 kph
- Range: 40-50 km
- Turn radius: 161 cm
- Ground clearance: 13cm
- Battery 12V/60-85A/H x 2pcs
- Front/Rear Wheel: 12" pneumatic

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au