



MERITS AURORA 2 HILL CLIMBER 4 WHEEL SCOOTER



SKU: S846C-HC-GR-B-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **Merits Aurora 2 Hill Climber Scooter** is purpose-built for users in need of powerful, reliable mobility—especially across inclines and uneven terrain. Its high-torque 4-pole motor and robust controller deliver exceptional hill-climbing performance without sacrificing agility. The 3-wheel design offers a tight turning radius, perfect for navigating indoors and in tight public spaces, while loop handlebars and intuitive controls make manoeuvring effortless. Enhanced stability and safety features—like the half-speed safety switch, emergency braking system, and adjustable suspension—ensure confidence and comfort during both everyday errands and outdoor excursions.

Features

- High-torque motor and upgraded controller provide superior hill-climbing capability
- Loop handlebars with finger or thumb throttle for ergonomic control
- Adjustable suspension system offers enhanced ride comfort and control
- Large 12" pneumatic tyres increase ground clearance for smoother travel
- Emergency braking system for added stopping security
- Half-speed safety switch automatically engages when turning to improve stability

- Convenient 12 V socket for charging devices or accessories

Specifications

- **Dimensions (4-wheel equivalent):** Approx. 146 cm (L) × 65 cm (W) × 123 cm (H)
- **Battery:** 12V/75A/H x 2pcs Ah
- **Turning Radius:** Approx. 161 cm
- **Ground Clearance:** ~13 cm
- **Range:** 40-50 km per charge
- **Max Speed:** 10 kph (limited to comply with Australian standards)
- **User Weight Capacity:** Up to 182 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au