



## **MERITS AURORA 2 HILL CLIMBER SCOOTER**

**\$5,999.23**



**SKU: S846C-HC-BR-H**

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits of Merits Aurora 2 Hill Climber Scooter**

The **Merits Aurora 2 Hill Climber Scooter** is a powerful and compact 3-wheel mobility scooter designed for users who require excellent hill-climbing capability and superior manoeuvrability. Ideal for both indoor use and outdoor travel on sloped terrain, the Aurora 2 provides smooth, stable performance with its robust motor and smart design.

Its 3-wheel configuration allows for a tighter turning radius, making it easy to navigate in shopping centres and narrow pathways. A padded high-back seat, intuitive controls, and front suspension ensure a comfortable ride, while the easy disassembly system makes it a practical choice for everyday mobility.

### **Key Features**

- High-torque motor for enhanced hill-climbing performance
- 3-wheel design for tight turning and easy navigation in confined spaces
- Padded swivel seat with adjustable armrests for comfort and access
- Front suspension for a smoother ride over uneven surfaces

- Disassembles easily for transport and storage
- Ideal for shopping, local errands, and hilly outdoor environments
- Simple, user-friendly tiller controls
- Front basket for carrying personal items

## **Specifications**

- Weight capacity: 182 kg
- Max speed: 12 kph (limited in Au to 10 km/h)
- Max speed (Hill Climber): 10 kph
- Range: 40-50 km
- Turn radius: 161 cm
- Ground clearance: 13cm
- Battery 12V/60-85A/H x 2pcs
- Front/Rear Wheel: 12" pneumatic