



MERITS AURORA S2 3-WHEEL SCOOTER - BRIGHT RED

\$5,999.99



SKU: S836C-BR

PRODUCT DESCRIPTION

Purpose & Benefits of Merits Aurora S2 3-Wheel Scooter - Bright Red

The **Merits Aurora S2 3-Wheel Scooter - Bright Red** is a compact, manoeuvrable mobility scooter designed for users who require comfort, control, and convenience in their everyday travels. Ideal for both indoor and outdoor use, this scooter's 3-wheel design offers a tight turning radius, making it perfect for navigating through shopping centres, narrow walkways, and crowded spaces.

Equipped with a high-back adjustable seat, easy-to-use controls, and a smooth suspension system, the Aurora S2 provides a stable and enjoyable ride. Its bright red finish adds a stylish touch, while the powerful motor and long-range battery ensure reliability for longer outings or daily errands.

Key Features:

- 3-wheel design for superior manoeuvrability in tight spaces
- Adjustable high-back seat for enhanced comfort and support
- Bright red colour for added visibility and style
- Suspension system for a smooth and stable ride
- Easy-to-use controls and delta tiller for intuitive steering

- Long-lasting battery with reliable range for daily use
- Suitable for indoor and outdoor environments
- Front basket for convenient personal storage

Specifications:

- Weight capacity: 182 kg
- Max speed: 12 kph (limited in Au to 10 km/h)
- Range: 40-50 km
- Turn radius: 161 cm
- Ground clearance: 13 cm
- Battery 12V/60-85A/H x 2pcs
- Front/Rear Wheel: 12" pneumatic
- Overall Dimension (3-Wheel) : 141 cm x 65 cm x 123 cm
- Seat Dimension : 46 cm x 40 cm x 51 cm