



MERITS P113 FOLD & GO COMPACT POWER WHEELCHAIR

\$3,142.80 - \$3,328.20 Price range: \$3,142.80 through \$3,328.20



SKU: P113

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Merits P113 Fold & Go compact power wheelchair folds in seconds for easy transportation or storage.

- Folds in one easy step
- Compact size means great maneuverability in tight spaces
- Removable battery box with battery charge
- indicator & extra charge port for charging off board
- Adjustable height footplate
- Adjustable height backrest
- Intelligent, regenerative, electromagnetic braking

SPECIFICATIONS

- Weight Capacity 114kg / 250lbs
- Max. Speed 6.4kph / 4mph
- Range 11km / 7mi. ; 18km / 11.25mi. (optional)
- Turn Radius 83cm / 33"
- Ground Clearance 3cm / 1.2"
- Motor DC 24V, 70W(Rating) / 210W(Max)
- Controller PG VR2 50A
- Battery 24V/ 11.4AH Li-Ion x 1pc
- Charger 2A off-board

- Gradient 6°
- Caster 6"(150×30) PU tire
- Drive Wheel 8"(200×50) solid tire
- Brake Intelligent, regenerative, electromagnetic brakes