




MERITS EXPLORER SCOOTER

~~\$6,145.00~~ Original price was:

~~\$6,145.00~~ \$5,530.50 Current price is: \$5,530.50.



SKU: S941A-BR-H

PRODUCT DESCRIPTION

Purpose & Benefits of Merits Explorer Scooter

The **Merits Explorer Scooter** is a heavy-duty 4-wheel mobility scooter designed for users who need reliable performance and comfort across a wide range of outdoor terrains. With a powerful motor, high ground clearance, and full suspension, the Explorer is built for extended travel and everyday independence in suburban and rural settings.

This scooter offers a spacious high-back seat with adjustable armrests, ergonomic controls, and a large front basket for convenience. Its robust construction and superior stability make it an excellent choice for users who require a durable and capable mobility solution for longer journeys and uneven surfaces.

Key Features

- Heavy-duty 4-wheel design for maximum stability
- Full suspension and large wheels for smooth travel over rough terrain
- High-back padded seat with swivel and adjustable armrests
- High ground clearance for outdoor accessibility
- Powerful motor and long-range battery for all-day travel
- LED lighting system for enhanced visibility and safety

- Ergonomic tiller with easy-to-use controls
- Ideal for independent outdoor mobility and extended travel

Specifications

- Max Speed 10 – 12 kph (we limit to 10 kph in Au)
- Range 40-50 km
- Turn Radius 167 cm
- Ground clearance 9 cm
- Motor DC 24V, 400 W (Rating)/1750 W (Max)
- Controller Dynamic Rhino2-12 0A
- Battery 12 V/44-75 A/H x 2 pcs
- Charger 5 A off-board
- Gradient 10°
- Front Wheel 12" pneumatic, Rear Wheel 12" pneumatic
- Brake Electromagnetic brake
- Overall Dimension : L x W x H 130 x 69 x 125 cm/51 x 27 x 49"
- Seat Dimension : W x D x H 46 x 40 x 51 cm/18 x 16 x 20"
- Seat Height from Ground : 63.5 cm/25"
- Seat Height from Deck : 43-48 cm/17-19"
- Wheelbase : 97 cm/38"