



MERITS FENDE MID-SIZE 4-WHEEL ELECTRIC SCOOTER



SKU: S946BR-H

PRODUCT DESCRIPTION

Purpose & Benefits of the Merits Fende Mid-Size 4-Wheel Electric Scooter

The Merits Fende Mid-Size 4-Wheel Electric Scooter is designed to provide individuals with a reliable and comfortable means of indoor and outdoor mobility. Its modern design, combined with advanced features, ensures users can navigate various terrains with ease and confidence.

Features

- **Stylish Design** – Features a modern aesthetic with a modesty panel, offering a smart and attractive appearance.
- **Adjustable Suspension** – Equipped with car-like front and rear suspension systems, delivering a smooth ride over varied terrains.
- **Comfortable Seating** – Includes a leather-feel, contoured high-back seat with slide, swivel, and recline functions, along with adjustable armrests for personalised comfort.
- **Enhanced Storage** – Unique battery storage design allows for greater floor space and provides a secure, lockable storage area for personal items.
- **LCD Dashboard** – Modern LCD offers clear and easy-to-read information, including battery gauge, speed, time, and temperature.
- **LED Lighting System** – Delivers brighter illumination with less power consumption, enhancing visibility and safety.
- **Safety Features** – Built-in speed reduction switch activates when turning, offering exceptional safety.

- **Convenient Power Access** – Includes an easy-to-access 12-volt socket suitable for charging devices like phones or GPS units.

Specifications

- **Weight Capacity:** Supports up to 136 kg.
- **Maximum Speed:** Up to 10 km/h.
- **Driving Range:** Approximately 30 km per charge.
- **Overall Dimensions:**
 - Length: 1200 mm
 - Width: 630 mm
 - Height: 1320 mm
- **Seat Dimensions:**
 - Width: 460 mm
 - Depth: 410 mm
 - Backrest Height: 460 mm
- **Turning Radius:** 1430 mm.
- **Ground Clearance:** 100 mm.
- **Climbing Grade:** Capable of handling inclines up to 12 degrees.
- **Total Weight:** 114.5 kg (including batteries).
- **Motor:** DC 24V, 295W (rated)/1150W (max).
- **Battery:** Four 12V/22Ah batteries.
- **Charger:** 5A off-board.
- **Wheels:** Both front and rear wheels are 10 inches (ca. 25 cm) in diameter.
- **Brakes:** Intelligent, regenerative, electromagnetic brakes.

The Merits Fende Mid-Size 4-Wheel Electric Scooter is an excellent choice for those seeking a balance of style, comfort, and functionality in their mobility solutions.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au