



# MERITS FENDE MID-SIZE 4-WHEEL ELECTRIC SCOOTER

\$5,282.00 Original price was: \$5,282.00.<u>\$4,754.00</u>Current price is: \$4,754.00.



**SKU:** S946BR-H

#### **PRODUCT DESCRIPTION**

## Purpose & Benefits of the Merits Fende Mid-Size 4-Wheel Electric Scooter

The Merits Fende Mid-Size 4-Wheel Electric Scooter is designed to provide individuals with a reliable and comfortable means of indoor and outdoor mobility. Its modern design, combined with advanced features, ensures users can navigate various terrains with ease and confidence.

#### **Features**

- **Stylish Design** Features a modern aesthetic with a modesty panel, offering a smart and attractive appearance.
- **Adjustable Suspension** Equipped with car-like front and rear suspension systems, delivering a smooth ride over varied terrains.
- **Comfortable Seating** Includes a leather-feel, contoured high-back seat with slide, swivel, and recline functions, along with adjustable armrests for personalised comfort.
- **Enhanced Storage** Unique battery storage design allows for greater floor space and provides a secure, lockable storage area for personal items.
- **LCD Dashboard** Modern LCD offers clear and easy-to-read information, including battery gauge, speed, time, and temperature.
- **LED Lighting System** Delivers brighter illumination with less power consumption, enhancing visibility and safety.
- **Safety Features** Built-in speed reduction switch activates when turning, offering exceptional safety.



 Convenient Power Access – Includes an easy-to-access 12-volt socket suitable for charging devices like phones or GPS units.

### **Specifications**

• Weight Capacity: Supports up to 136 kg.

• Maximum Speed: Up to 10 km/h.

• **Driving Range:** Approximately 30 km per charge.

• Overall Dimensions:

Length: 1200 mmWidth: 630 mmHeight: 1320 mm

• Seat Dimensions:

Width: 460 mmDepth: 410 mm

Backrest Height: 460 mm
Turning Radius: 1430 mm.
Ground Clearance: 100 mm.

• Climbing Grade: Capable of handling inclines up to 12 degrees.

Total Weight: 114.5 kg (including batteries).
Motor: DC 24V, 295W (rated)/1150W (max).

• Battery: Four 12V/22Ah batteries.

• Charger: 5A off-board.

• Wheels: Both front and rear wheels are 10 inches (ca. 25 cm) in diameter.

• **Brakes:** Intelligent, regenerative, electromagnetic brakes.

The Merits Fende Mid-Size 4-Wheel Electric Scooter is an excellent choice for those seeking a balance of style, comfort, and functionality in their mobility solutions.