



MERITS HEAVY DUTY ROLLATOR

\$319.00



SKU: RC-SC5089

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Merits Heavy Duty Rollator is a durable and reliable rollator that is designed to support users up to 136kg. It features a sturdy aluminium frame and a variety of features that make it comfortable and easy to use.

FEATURES

- Sturdy aluminium frame
- Adjustable height handlebars
- Padded seat and backrest
- Basket for carrying items
- Loop brakes for slowing down and stopping
- 8" puncture-proof polyurethane wheels
- Swivel front casters for maneuverability
- Easy to fold and transport

The Merits Heavy Duty Rollator is ideal for people of all ages and abilities who need a durable and reliable rollator. It is especially well-suited for people who have difficulty with balance or who need a rollator to support their weight.

SPECIFICATIONS

- Depth 700mm
- Height - 800 mm to 915mm
- Load Capacity - 182kg

- Seat Depth – 335mm
- Seat Height – 525mm
- Seat width – 465mm
- Weight – 9kg
- Width – 660mm