



## **MERITS HEAVY DUTY ROLLATOR**



**SKU:** RC-SC5089

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The Merits Heavy Duty Rollator is a durable and reliable rollator that is designed to support users up to 136kg. It features a sturdy aluminium frame and a variety of features that make it comfortable and easy to use.

### **FEATURES**

- Sturdy aluminium frame
- Adjustable height handlebars
- Padded seat and backrest
- Basket for carrying items
- Loop brakes for slowing down and stopping
- 8" puncture-proof polyurethane wheels
- Swivel front casters for maneuverability
- Easy to fold and transport

The Merits Heavy Duty Rollator is ideal for people of all ages and abilities who need a durable and reliable rollator. It is especially well-suited for people who have difficulty with balance or who need a rollator to support their weight.

### **SPECIFICATIONS**

- Depth 700mm
- Height - 800 mm to 915mm
- Load Capacity - 182kg

- Seat Depth – 335mm
- Seat Height – 525mm
- Seat width – 465mm
- Weight – 9kg
- Width – 660mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)