




MERITS INTERCEPTOR S2 4 WHEEL SCOOTER

\$8,298.71



SKU: S940A-BR-H

PRODUCT DESCRIPTION

Purpose & Benefits of Merits Interceptor S2 4 Wheel Scooter

The **Merits Interceptor S2 4 Wheel Scooter** is a high-performance mobility scooter built for users who value speed, stability, and advanced comfort features. Designed for long-distance travel and outdoor use, the Interceptor S2 delivers exceptional performance with a powerful motor, full suspension system, and large pneumatic tyres for a smooth ride on varied terrain.

Ideal for active individuals seeking greater independence, this scooter features a high-back captain's seat, ergonomic tiller controls, and a modern digital dashboard. With its strong build and stylish design, the Interceptor S2 combines functionality and comfort for daily errands, recreational outings, and extended trips.

Key Features

- Heavy-duty 4-wheel design for enhanced stability and safety
- Full suspension and large pneumatic tyres for a smooth outdoor ride
- High-back captain's seat with adjustable headrest and armrests
- Digital dashboard with speed, battery level, and trip information
- Powerful motor suitable for longer distances and varied terrain
- Front and rear LED lighting for increased visibility and safety

- Ergonomic tiller with easy-grip controls
- Ideal for users who require a premium, road-ready scooter

Specifications

- Weight capacity: 204 kg
- Max speed: 10-15 km/h (Limited to 10 km/h in Au)
- Range: 40 – 50 km
- Total Scooter Weight: 148 kg
- Turning radius: 180 cm
- Overall length: 154 cm, Width: 72 cm, Height: 145.5 cm
- Seat width: 20" (51 cm)
- Colours: Red, Silver, Pearl Gold