



#### **MERITS INTERCEPTOR SCOOTER - BRIGHT RED**

<del>\$7,956.30</del> Original price was: \$7,956.30.<u>\$7,161.00</u>Current price is: \$7,161.00.



**SKU:** S940A-1BR

#### PRODUCT DESCRIPTION

### Purpose & Benefits of Merits Interceptor Scooter - Bright Red

The **Merits Interceptor Scooter - Bright Red** is a high-performance mobility scooter designed to provide comfort, safety, and style for individuals who require reliable, long-distance travel support. With its sporty design and robust build, the Interceptor is ideal for active users seeking both functionality and aesthetics.

Equipped with a full suspension system, high-back captain's seat, and advanced digital controls, the Interceptor ensures a smooth and enjoyable ride over various terrains. Its powerful motor and long-range battery make it perfect for daily errands, outings, and independent travel.

Designed with user comfort and convenience in mind, the Interceptor offers ample legroom, adjustable features, and a modern LED lighting system for enhanced visibility and safety.

## **Key Features:**

- Sporty design with bright red finish for enhanced visibility and style
- Full suspension and pneumatic tyres for a smooth ride
- High-back adjustable captain's seat for comfort and support
- Digital dashboard display with speed, battery, and trip info



- Powerful motor with a long-range battery for extended travel
- Front and rear LED lights and indicators for added safety
- Large rearview mirrors and ergonomic handlebars
- Suitable for outdoor use on varied terrain

# **Specifications:**

Overall Length	1540 mm
Overall Width	720 mm
Overall Height	1450 mm
Wheelbase Length	1130 mm
Seat Width	510 mm
Seat Depth	460 mm
Min Floor To Seat Height	430 mm
Max Floor To Seat Height	500 mm
Min Seat Height from Deck	600 mm
Max Seat Height from Deck	680 mm
Turning Radius	1800 mm
Backrest Height	630 mm
min-range KM	40 km
max-range KM	60 km
Max Speed	10 km/h
Front Whool	1/1" ppour

Front Wheel 14" pneumatic Rear Wheel 16" pneumatic

Gradient 10°

Ground Clearance 130 mm

Battery 12V/75A/H x 2pcs Ah
Motor Size 650W/2200W hp
Controller PG S-Drive 200A

Charger 8Amp A
Product Weight 152 kg
Total Weight Without Battery 102 kg
Max User Weight 204 kg



Brand Merits