




MERITS INTERCEPTOR SCOOTER - BRIGHT RED



SKU: S940A-1BR

PRODUCT DESCRIPTION

Purpose & Benefits of Merits Interceptor Scooter - Bright Red

The **Merits Interceptor Scooter - Bright Red** is a high-performance mobility scooter designed to provide comfort, safety, and style for individuals who require reliable, long-distance travel support. With its sporty design and robust build, the Interceptor is ideal for active users seeking both functionality and aesthetics.

Equipped with a full suspension system, high-back captain's seat, and advanced digital controls, the Interceptor ensures a smooth and enjoyable ride over various terrains. Its powerful motor and long-range battery make it perfect for daily errands, outings, and independent travel.

Designed with user comfort and convenience in mind, the Interceptor offers ample legroom, adjustable features, and a modern LED lighting system for enhanced visibility and safety.

Key Features:

- Sporty design with bright red finish for enhanced visibility and style
- Full suspension and pneumatic tyres for a smooth ride
- High-back adjustable captain's seat for comfort and support
- Digital dashboard display with speed, battery, and trip info

- Powerful motor with a long-range battery for extended travel
- Front and rear LED lights and indicators for added safety
- Large rearview mirrors and ergonomic handlebars
- Suitable for outdoor use on varied terrain

Specifications:

Overall Length	1540 mm
Overall Width	720 mm
Overall Height	1450 mm
Wheelbase Length	1130 mm
Seat Width	510 mm
Seat Depth	460 mm
Min Floor To Seat Height	430 mm
Max Floor To Seat Height	500 mm
Min Seat Height from Deck	600 mm
Max Seat Height from Deck	680 mm
Turning Radius	1800 mm
Backrest Height	630 mm
min-range KM	40 km
max-range KM	60 km
Max Speed	10 km/h
Front Wheel	14" pneumatic
Rear Wheel	16" pneumatic
Gradient	10°
Ground Clearance	130 mm
Battery	12V/75A/H x 2pcs Ah
Motor Size	650W/2200W hp
Controller	PG S-Drive 200A
Charger	8Amp A
Product Weight	152 kg
Total Weight Without Battery	102 kg
Max User Weight	204 kg

Brand

Merits

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au