



MERITS INTERCEPTOR SCOOTER



SKU: N/A

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

This large sporty streamline design scooter has full suspension with shock absorbers for a more comfortable ride. Your battery gauge, temperature, accumulated usage hours and maintenance messages are all displayed on the LCD digital dashboard. This scooter also has a handy storage space at the tiller and a comfortable seating system.

FEATURES

- Brand new streamline design with a sporty feel
- Fully adjustable suspension system that delivers a smooth and comfortable ride
- Modern LED lighting system delivers brilliant light with less power consumption
- Digital LCD dashboard display including a time display, battery gauge, temperature, maintenance messages, speed display and odometer
- Extra storage space beneath tiller
- Built-in half speed switch which operates automatically when turning to provide extra safety
- Merits 4PT (4 pole technology) and large controller, delivering more power for climbing steep inclines
- New easy to access 12 volt socket

SPECIFICATIONS

Overall Length: 154 cmOverall Width: 72 cmOverall Height: 145.5 cm

• Seat height from ground: 60-68 cm



• Seat height from deck: 43-50 cm

• Wheel base: 113 cm

• Total Capacity: 204 kg (450 lbs)

• Max Speed: 10-15 km

• Max Range per charge: 40-60 km

• Turn Radius: 180 cm

• Ground clearance: 13 cm

Motor: DC 24 V, 650W (rating) / 220 W (max)

Controller: PG S-drive 200ABatteries: 12V , 60-100AH x 2

Charger: 8A off-boardGradient: 10 degrees

• Tires: Pneumatic

• front 14"

∘ rear 16"

Brakes: Intelligent, regenerative, electromagnetic brakes

• Total Weight with Batteries: 148 kg (325 lbs)

Buy This Merits Interceptor Scooter online at MediMart using Afterpay or <u>NDIS method</u> today! Expert <u>mobility equipment suppliers</u> in Australia.

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

<u>brighton@medimart.com.au</u>

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 319903 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au



NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au